

Product Information : StemFit Basic03

1. Introduction

StemFit Basic03 is an animal-origin free, defined medium for human pluripotent stem cell culture. This medium is capable of reprogramming, maintaining and differentiating human ES and iPS cells under feeder-free culture conditions.

2. Materials Provided 1





Component	Volume	Storage
Liquid A	400ml	Store at below -20 °C
Liquid B	100ml	Store at below -20 °C



3. Media Preparation

StemFit Basic03 is provided frozen as a 2 component set containing “Liquid A” and “Liquid B” and can be stored at below -20 °C until use. Use sterile techniques to prepare StemFit Basic03 medium.

- 1) Before use, thaw the frozen “Liquid A” and “Liquid B” with occasionally mixing at room temperature (15-25 °C). 2
- 2) Aseptically mix medium components by adding the full volume of “Liquid B” to “Liquid A”. Mix thoroughly. If prepared aseptically, StemFit Basic03 is ready for use.
- 3) Upon thawing, StemFit Basic03 medium may be aseptically aliquoted and stored at below -20 °C. 3 4
- 4) Before use, warm aliquots to room temperature and use immediately.

<div style="background-color: #f4a460; padding: 10px; text-align: center;"> 1  bFGF not contained </div> <div style="padding: 10px;"> <p>This product does not contain basic fibroblast growth factor (bFGF)</p> </div>	<div style="background-color: #f4a460; padding: 10px; text-align: center;"> 2  × 37 °C </div> <div style="padding: 10px;"> <p>Do not thaw, “Liquid B” at 37 °C, as it accelerates the degradation of medium ingredients.</p> </div>
<div style="background-color: #f4a460; padding: 10px; text-align: center;"> 3  2 weeks </div> <div style="padding: 10px;"> <p>Thawed StemFit Basic03 medium may be stored at 2-8 °C for up to two weeks.</p> </div>	<div style="background-color: #f4a460; padding: 10px; text-align: center;"> 4  × Light </div> <div style="padding: 10px;"> <p>We recommend storing the medium protected from light</p> </div>

4. Precaution and Disclaimer

StemFit Basic03 is for research use only and is not for diagnostic use, therapeutic use (including clinical research), or as a food.

5. Contact the following department for product information:

Amino Acids Dept., AminoScience Division, AJINOMOTO CO., INC.

1-15-1 Kyobashi, Chuo-ku, Tokyo 104-8315, Japan

TEL : +81-3-5250-5040 E-mail: stemfit@ajinomoto.com

For Research Use Only. Not Intended For Human or Animal Diagnostic or Therapeutic Uses