FabRICATOR®Z

FOR RESEARCH USE ONLY

www.genovis.com

-20°C



SmartEnzymes™



FabRICATOR®Z

INSTRUCTIONS FOR PRODUCT

FabRICATOR®Z 2000 units (A0-FRZ-020) Digestion of up to 2 mg mouse IgG2a or mouse IgG3

QUICK GUIDE

Prepare FabRICATOR®Z

Reconstitute FabRICATOR Z in $50 \ \mu l \ ddH_2O$ to a concentration of $40 \ units/\mu l$.



2 Add FabRICATOR®Z

Add 1 unit FabRICATOR Z / 1 µg IgG1.



3 Digestion
Incubate for 2 h² at 37°C.



PRODUCT DESCRIPTION

FabRICATOR Z is an enzyme for preparation of F(ab')2 from mouse subclasses IgG2a and IgG3. FabRICATOR Z digests mouse IgG2a and IgG3 at one specific site below the hinge region, resulting in precise F(ab')2 and Fc fragments (Table 1). Most mouse IgG2a and IgG3 are digested by FabRICATOR Z within 2 hours, although the enzyme:Ab ratio and/or digestion time may need to be optimized for certain antibodies. For digestion of human IgGs, FabRICATOR® is recommended.

FabRICATOR Z digests IgG at one specific site and there is no risk of overdigestion if the incubation time is prolonged. FabRICATOR Z digests IgG at mild conditions, thus preserving the immunoreactivity. FabRICATOR Z digests IgG in sodium phosphate buffer with pH ranging from 6.0-8.0, but is most active at pH 6.5 and at low NaCl concentrations.

FabRICATOR Z is derived from S. equi ssp. zooepidemicus and expressed in E. coli. The enzyme contains a His-tag and the molecular weight is 36 442 Da.

Immobilized FabRICATOR Z is available in convenient spin column format, FragIT™Z, for generation of fragments without FabRICATOR Z in the final sample preparation.

FabRICATOR® Z

Table 1. FabRICATOR Z digestion sites.

Mouse IgG subclass	FabRICATOR Z digestion site
lgG2a	CPAPNLLG / GPSVF
lgG3	CPPGNILG / GPSVF

Unit Definition

One unit FabRICATOR Z digests \geq 95% of 1 μg mouse IgG2a when incubated in 10 mM sodium phosphate, 10 mM NaCl, pH 6.5 at 37°C for 2 h.

Content and Storage

FabRICATOR Z is supplied lyophilized in 50 mM Tris, 150 mM NaCl, pH 7.6, with no preservatives added.

FabRICATOR Z is shipped at ambient temperature and should be stored at -20°C upon arrival.

After reconstitution, FabRICATOR Z is stable for 1 month at +4-8°C.

FabRICATOR Z is for R&D use only.

DETAILED PROTOCOL

Additional Materials Required

Digestion buffer³: 10 mM sodium phosphate, 10 mM NaCl, pH 6.5.

Preparation of IgG

Prepare the IgG in digestion buffer³. The final IgG concentration in the reaction should be 0.5-10 mg/ml.

FabRICATOR®Z

Antibody Subunit Generation

- 1 Prepare FabRICATOR®Z

 Reconstitute FabRICATOR Z in 50 μI

 ddH₂O to a concentration of 40 units/μI.
- 2 Add FabRICATOR®Z

 Add 1 unit FabRICATORZ / 1 μg IgG¹.

 The final IgG concentration should be 0.5-10 mg/ml.
- 3 Digestion
 Incubate for 2 h² at 37°C

Notes

- The enzyme:antibody ratio may need to be increased for individual antibodies.
- The incubation time may need to be prolonged for individual antibodies.
- pH and NaCl concentration can be increased but this may require increased enzyme concentration and/or prolonged incubation time

Fabricator®Z

Quality Control

FabRICATOR Z is tested to meet the specifications and lot-to-lot consistency.

FabRICATOR Z is tested for absence of microbial contamination with blood agar plates, Sabouraud dextrose agar plates and fluid thioglycollate medium.

Related Products

FragIT™ Z

Immobilized FabRICATOR Z, digestion of mlgG2a and mlgG3

FragIT™ Z kit

Digestion of mlgG2a or mlgG3 and purification of F(ab')2 and Fc fragments

FabRICATOR® Z Limited Use Label License: Research Use Only

All rights reserved. The trademark FabRICATOR®Z is the property of Genovis AB. FabRICATOR® is duly registered in Australia, China, EU, Japan, South Korea, USA and Canada.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

All goods and services are sold subject to Genovis' General Terms and Conditions of Sale.

©2017 Genovis AB.


OTHER PRODUCTS



Digestion of IgG and Purification of F(ab')2 and Fc Fragments

FragIT kit consists of an IgG digestion column, FragIT, and an affinity purification column, CaptureSelect™*. FragIT is a resin with FabRICATOR® enzyme covalently coupled to agarose beads for subunit fragmentation of IgG to generate F(ab')2 and Fc fragments. After digestion, the fragments can be purified using the CaptureSelect column supplied in the kit.

- Digestion of IgG on a column.
- · Purification of F(ab')2 and Fc fragments.




*Made with Thermo Scientific™ CaptureSelect™ resin from Thermo Fisher Scientific Inc. and its subsidiaries. Thermo Scientific and CaptureSelect are trademarks of Thermo Fisher Scientific Inc. and its subsidiaries.

FabULOUS®

Digestion of IgG from Many Different Species and Subclasses

FabULOUS is a cysteine protease that digests IgG in the upper hinge region, generating intact Fab and Fc fragments.

- Active on several IgG species, including human, mouse, rat and goat.
- Rapid a 60 minute protocol.
- FabULOUS Fab kit is available for rapid and easy purification of Fab fragments from mouse IgG.





US & Canada

Genovis Inc. 245 First Street, Suite 1800 Cambridge, MA 02142 USA

Customer service: 001 (617)-444-8421 Order phone (toll free): 001 (855)-782-0084

Order fax: 001 (858)-524-3006 Email: orders.us@genovis.com

EMEA & Asia

Genovis AB Box 790 SE-220 07 Lund Sweden

Customer service: 0046 (0)46 10 12 30 Order phone: 0046 (0)46 10 12 30 Order fax: 0046 (0)46 12 80 20 Email: order@genovis.com