



FucosEXO™




Hydrolysis of α1-2,3,4 linked Fucose

SmartEnzymes™







FucosEXO™



Hydrolysis of α 1-2,3,4 linked fucose on glycoproteins or oligosaccharides.

FucosEXO is a mix of α -fucosidases for efficient hydrolysis of α 1-2, α 1-3 and α 1-4 linked fucose residues on native *N*- and *O*-glycosylated proteins or free oligosaccharides. The two enzymes are modified with His-tags and expressed in *E. coli*. The molecular weight of the enzymes are 87 kDa and 64 kDa, respectively. FucosEXO is active in neutral pH and does not require any co-factors or special buffers. The exoglycosidase activity of FucosEXO allows for the complete hydrolysis of α 1-2,3,4 linked fucose on up to 2 mg of glycoprotein. FucosEXO is a valuable tool for glycan structure analysis on *N*- and *O*-glycoproteins or in exoglycosidase arrays on oligosaccharides.

-  Hydrolyzes α 1-2,3,4 linked Fucose on native glycoproteins
-  α 1-2,3,4 linked Fucose
-  1-2 h incubation
-  No co-factors required


Key Characteristics

- ▶ Efficient hydrolysis of α 1-2,3,4 linked fucose
- ▶ For native glycoproteins or oligosaccharides
- ▶ Activity without the need for co-factors

Specific Hydrolysis of Fucose Linkages

Fucose can be attached to both *N*- and *O*-glycans. The α 1-2, α 1-3 and α 1-4 linkages most commonly occur in *O*-glycans or as antenna fucosylation of *N*-glycans, whereas α 1-6 linked fucose is found as a modification of the *N*-glycan core. We measured the release of fucose using FucosEXO from a panel of oligosaccharide substrates representing different linkages. FucosEXO efficiently releases α 1-2, α 1-3 and α 1-4 linked fucose residues (Fig. 1).


Figure 1. The substrate specificity of FucosEXO was analyzed on equal molar amounts of synthetic oligosaccharides incubated with FucosEXO for 30 min at 37°C and the amount of released fucose was measured spectrometrically using a L-Fucose Assay Kit (Megazymes).



Defucosylation of Native Glycoproteins

Fucosylation of O-glycans is involved in synthesis of functionally important glycan epitopes such as blood group antigens and the Lewis structures. The analysis of glycoproteins modified with such complex glycans can be challenging and requires specific and efficient enzymatic tools. We tested FucosEXO on a glycoengineered TNFR protein carrying up to 11 O-glycans decorated with α 1-2 and α 1-3 linked fucose and compared the activity to a commercially available α -fucosidase. Within 1 hour, complete removal of fucose on the TNFR protein was achieved with FucosEXO, while treatment with the other fucosidase only had a minor impact on fucoses (Fig. 2). Hydrolysis of all fucose residues from the TNFR protein was also achieved using Immobilized FucosEXO.

Figure 2. Glycoengineered TNFR with core 1 and core 2 O-glycans decorated with both 1-2 and 1-3 fucose was incubated with bovine kidney fucosidase (1h, 37°C), FucosEXO (1h, 37°C) or Immobilized FucosEXO (1h, RT). The resulting protein was analyzed by reverse phase LC-MS on a Waters BioAccord system equipped with a Waters BioResolve RP mAb column (2.1 x 50 mm).



No Additives or Co-factors Required

FucosEXO is active in physiological buffers at a neutral pH and the activity is not dependent on any co-factors, making FucosEXO compatible with most samples without the risk of impairing LC-MS analysis (Table 1). Depending on the nature of the substrate, FucosEXO hydrolyzes terminal α 1-2,3,4 linked fucose residues on glycoproteins within 1 to 2 hours, while longer incubation may be required for very complex samples. FucosEXO is also available immobilized on agarose beads in a spin column format for convenient processing without residual enzyme in the final preparation.

Table 1. Key characteristics of FucosEXO.

Enzyme Feature	FucosEXO
Incubation time	1-2 hours
pH range	6.0 - 8.0
MS compatibility	Yes
Special buffers	No
Co-factors	No
Additives or BSA	No

FucosEXO™



Product ID	Description	EUR	USD
G1-FM1-020	FucosEXO, 2000 units	750	845
G1-FM6-025	Immobilized FucosEXO, 5 x 0.5 mg	795	995
G1-FM6-050	Immobilized FucosEXO, 10 x 0.5 mg	1,295	1,745

Legal and Disclaimers

All rights reserved. Genovis products may be covered by one or more patents, trademarks and copyrights owned or controlled by Genovis AB.
For more information about commercial rights, please contact the Genovis team at licensing@genovis.com.
Genovis products are intended for research use only. They are not intended to be used for therapeutic or diagnostic purposes in humans or animals.

©2021 Genovis AB

US & Canada

Genovis Inc.
245 First Street, Suite 1800
Cambridge, MA 02142
USA

Customer service: 001 (617)-444-8421
Order phone (toll free): 001 (855)-782-0084
Order fax: 001 (858)-524-3006
Email: orders.us@genovis.com

EMEA & Asia

Genovis AB
Box 790
SE-220 07 Lund
Sweden

Customer service: +46 46 10 12 30
Order phone: +46 46 10 12 30
Order fax: +46 46 12 80 20
Email: order@genovis.com



info@genovis.com | www.genovis.com