

# MSOT MULTISPECTRAL OPTOACOUSTIC TOMOGRAPHY ACUITY



## MSOT Acuity | Acuity Echo

- Clinical imaging system for research and diagnosis
- Tomographic handheld real-time imaging
- Anatomical, functional and molecular contrast

## MSOT TECHNOLOGY

### Imaging sequence:

- Illumination of tissue with laser pulses at multiple wavelengths
- Detection of induced ultrasound pressure waves
- Spectral unmixing to analyze individual absorbers

### Technology benefits:


Combines the molecular specificity of optical imaging with the depth and spatiotemporal resolution of ultrasound

#### • Molecular specificity:

Identify and quantify disease-related biomarkers, revealing endogenous absorbers and injected probes

#### • Depth & spatiotemporal resolution:

Acquire soft tissue images, with a spatial resolution of up to 80  $\mu\text{m}$ , *in vivo* and in real time



## i OPUS: hybrid OPTOacoustic & UltraSound imaging

Ultrasound imaging integrated in the MSOT Acuity Echo enables the visualization of both tomographic optoacoustic and ultrasound information at the same time, thus providing additional and complementary information on tissue morphology.

## MSOT VS. OTHER IMAGING MODALITIES

Imaging is an essential tool for medical diagnosis. Various technologies have evolved over time, each one providing specific benefits – but also limitations. Innovation in biomedical imaging facilitates new avenues for diagnosis and treatment of diseases.

**X-ray, ultrasound** and magnetic resonance imaging (**MRI**) provide anatomical information at high spatial resolution but with limited molecular specificity and sensitivity. On the other hand, positron emission tomography (**PET**) is a molecular imaging modality with high sensitivity, but suffers from low spatial resolution.


Modality	Spatial resolution	Temporal resolution	Sensitivity	Cost and infrastructure	Burden for patients
MSOT	●	●	○	●	●
X-ray	●	○	–	●	–
Ultrasound	●	●	–	●	●
MRI	●	○	–	–	○
PET/SPECT	–	–	●	–	–

● = favorable   ○ = medium   – = unfavorable

**MSOT** – Multispectral Optoacoustic Tomography is a novel imaging technology. Besides its high spatio-temporal resolution and sensitivity for optical contrast, MSOT imaging comes at relatively low cost.

Additionally, MSOT imaging poses no significant burden – particularly no ionizing radiation – on patients and users.

## SYSTEM COMPONENTS



### MSOT scanner hardware

- Mobile imaging platform with integrated high-performance workstation and adjustable screen
- Graphical user interface controlled via touchscreen and customized keyboard
- Detectors for different applications with variable center frequency, geometry, and size

### ViewMSOT™ software

- Intuitive software design
- Patient and study data administration
- Automated spectral unmixing and signal quantification
- Wide range of tools for data analysis and export
- Remote analysis support

## SYSTEM USE



### System design allows intuitive use:

- Workflow comparable to that of conventional ultrasound
- Handheld detectors for easy patient access
- Real-time spectral analysis of tissue chromophore distribution
- One click acquisition of single images and image sequences
- Detailed analysis and comparison of multiple images
- Export of image data and standard examination reports



## ACCESSORIES AND OPTIONS

### Standard accessories:

- Laser safety goggles
- User kit, incl. ultrasound gel, disinfectant wipes, sterile covers
- Calibration phantom

### Options:

- 2D and 3D detectors customized for particular application needs (2.5-10 MHz, 1-512 elements)
- Advanced software tools for image analysis
- Device for imaging of tissue specimens



## VISUALIZED CHROMOPHORES

MSOT detects and visualizes signals that represent the spectrally distinct absorbance of chromophores in tissue. Preclinical and clinical studies have proven the effective detection of endogenous chromophores such as hemoglobin, melanin and lipids as well as that of commonly used and clinically approved optical contrast agents such as ICG and methylene blue.


## CASE STUDIES

MSOT has already been applied in a variety of clinical applications, including the examples listed below:

### Nodal assessment in malignant melanoma<sup>1</sup>

MSOT has shown the potential for non-radioactive localization of sentinel lymph nodes (SLNs) by tracking the lymphatic drainage of ICG in melanoma and other cancers where SLN biopsy is a routine diagnostic procedure.


In melanoma, the absence of a melanin signal in the SLN can potentially rule out the presence of metastasis, and thereby spare many patients an unnecessary surgical procedure.



### Peripheral vascular diseases (PVD)<sup>2</sup>

MSOT has the ability to image vessels with a diameter of less than 100  $\mu\text{m}$  in real time and contrary to other imaging modalities such as X-ray, CT and MRI, without the need for contrast agents.


Images acquired by MSOT can show vascular malformation and blood/tissue oxygenation state, for example in PVD patients, thereby supporting the assessment of disease progression.



### Thyroid nodules<sup>3</sup>

Only 5% of all biopsied thyroid nodules are malignant. Optoacoustic tomography could help to distinguish malignant and benign nodules non-invasively by determining perfusion and tissue oxygenation of the nodules and showing vascularization patterns.


Similar methods are used to assess malignancy of other cancerous lesions, e.g. in breast cancer patients.



### Alopecia<sup>4</sup>

Monitoring of viable hair follicles over time is crucial for understanding hair pathophysiology and discovery of novel therapies for alopecia, but also other follicular diseases.

MSOT provides structural information and functional properties of the entire hair follicle, such as vascularization and oxygenation of the bulb and the lipid volume of the associated sebaceous glands.



1) Stoffels I et al., **Metastatic status of sentinel lymph nodes in melanoma determined noninvasively with multispectral optoacoustic imaging**, *Sci Transl Med*. 2015 Dec 9;7(317):317ra199. DOI: 10.1126/scitranslmed.aad1278.

2) Deán-Ben XL and Razansky D, **Functional optoacoustic human angiography with handheld video rate three dimensional scanner**, *Photoacoustics*. 2013 Nov 12;1(3-4):68-73. DOI: 10.1016/j.pacs.2013.10.002.

3) Dogra VS et al., **Preliminary results of ex vivo multispectral photoacoustic imaging in the management of thyroid cancer**, *AJR Am J Roentgenol*. 2014 Jun;202(6):W552-8. DOI: 10.2214/AJR.13.11433.


4) Ford SJ et al., **Structural and functional analysis of intact hair follicles and pilosebaceous units by volumetric multispectral optoacoustic tomography**, *J Invest Dermatol*. 2015 Dec 29. DOI: 10.1016/j.jid.2015.09.001.

# TECHNICAL SPECIFICATIONS

Specifications	MSOT Acuity <sup>1</sup>	MSOT Acuity Echo <sup>2</sup>	
<b>Image acquisition</b>			
Image rate (live display)	up to 50 fps		
Acquisition time	< 10 ms (single wavelength)   < 100 ms (multispectral)		
<b>Key system components</b>			
<b>SpectraPULSE™ illumination system</b>			
Laser wavelength spectrum	Standard: 680-980 nm   Optional: 670-1250 nm		
Pulse repetition rate	10-50 Hz (adjustable)		
Maximum pulse energy @ 750 nm	30 mJ		
PulseCTRL™	Pulse energy control, laser performance monitoring		
Wavelength tuning	Tuning time: < 10 ms   Minimum step size: 1 nm		
<b>RapidSCAN™ data acquisition electronics</b>			
Channels for simultaneous acquisition	up to 512		
Sampling rate	up to 40 MS/s		
<b>OPUS™ hybrid ultrasound mode</b>			
Transmit frequency	n/a	2-8 MHz	
Advanced imaging methods	n/a	Synthetic aperture beamforming, spatial compounding	
<b>ViewMSOT™</b>			
Data management	Patient/study data administration, creation of reports, data export		
Data acquisition	2D/3D, single images, image sequences, multispectral data sets		
Data processing	Image reconstruction, spectral unmixing, signal quantification		
Data analysis	Measurements, spectral analysis, color maps, image filters		
<b>General technical specifications</b>			
Hardware scanner console	Intel Core i7/Xeon, 32 GByte RAM, 4 TByte HDD data storage, 24" TFT touchscreen		
Operating system	Windows 8 embedded		
External Interfaces	1 GBit Ethernet, remote interlock connector		
Dimensions (width x depth x height   weight)	73 x 91 x 152 cm   290 kg		
Power	16A/230VAC, 50/60 Hz		
Laser classification	Class 4		
<b>TomoARC™ detectors</b>			
	<b>Standard 2D</b>	<b>Standard 3D</b>	<b>Optional detectors</b>
Angular coverage	125°	110°	90-180°
Center frequency	3 MHz	2.5 MHz	2.5-10 MHz   > 50%
Number of elements	256	384	1-512
Field of view	25 mm	15 x 15 mm	up to 30   20 x 20 mm
Depth penetration (MSOT)	up to 30 mm	up to 25 mm	up to 30 mm
Maximum resolution	275 µm	390 µm	80-390 µm

1) MSOT Acuity will be CE marked in 2016. Please contact us for further information.

2) MSOT Acuity Echo is currently only available as a research platform.



## SELECTION OF RELATED PUBLICATIONS

- Stoffels I et al.  
**Metastatic status of sentinel lymph nodes in melanoma determined noninvasively with multispectral optoacoustic imaging,**  
Sci Transl Med. 2015 Dec 9;7(317):317ra199.
- Taruttis A and Ntziachristos V.  
**Advances in real-time multispectral optoacoustic imaging and its applications,**  
Nat. Photonics. 2015;9:219–227.
- Zackrisson S et al.  
**Light in and sound out: emerging translational strategies for photoacoustic imaging,**  
Cancer Res. 2014 Feb 15;74(4):979-1004.
- Garcia-Urbe A et al.  
**Dual-Modality Photoacoustic and Ultrasound Imaging System for Noninvasive Sentinel Lymph Node Detection in Patients with Breast Cancer,**  
Sci Rep. 2015 Oct 29;5:15748.
- Dogra VS et al.  
**Preliminary results of ex vivo multispectral photoacoustic imaging in the management of thyroid cancer,**  
AJR Am J Roentgenol. 2014 Jun;202(6):W552-8.
- Buehler A et al.  
**Real-time handheld multispectral optoacoustic imaging,**  
Opt Lett. 2013 May 1;38(9):1404-6.
- Deán-Ben XL and Razansky D.  
**Functional optoacoustic human angiography with handheld video rate three dimensional scanner,**  
Photoacoustics. 2013 Nov 12;1(3-4):68-73.
- Razansky D et al.  
**Volumetric real-time multispectral optoacoustic tomography of biomarkers,**  
Nat Protoc. 2011 Jul 7;6(8):1121-9.
- Bayer CL et al.  
**Photoacoustic imaging: a potential tool to detect early indicators of metastasis,**  
Expert Rev Med Devices. 2013 Jan;10(1):125-34.
- Taruttis A et al.  
**Mesosopic and macroscopic optoacoustic imaging of cancer,**  
Cancer Res. 2015 Apr 15;75(8):1548-59.
- Mehrmohammadi M et al.  
**Photoacoustic Imaging for Cancer Detection and Staging,**  
Curr Mol Imaging. 2013 Mar;2(1):89-105.