

# OmniGLYZOR™

---

FOR RESEARCH  
USE ONLY  
[www.genovis.com](http://www.genovis.com)

STORE CONTENT  
AT DIFFERENT  
TEMPERATURES  
(See page 6)



## SmartEnzymes™

---



GENOVIS

---

## INSTRUCTIONS FOR PRODUCTS

---

**OmniGLYZOR, 5 × 50-100 µg** (G3-OM6-005)  
Deglycosylation of 5 × 50-100 µg glycoprotein

**OmniGLYZOR, 10 × 50-100 µg** (G3-OM6-010)  
Deglycosylation of 10 × 50-100 µg glycoprotein

### Quick Guide

- The Quick Guide (p.3) is intended for experienced users. First time users are recommended to follow the detailed protocol (p. 8).
- Use lids and bottom caps during the incubation.
- Before centrifugation, remove the bottom cap and slightly open the lid.
- If denaturing conditions are needed for PNGase F, an additional step is performed with lyophilized PNGase F in the presence of *RapiGest*™ SF Surfactant\*. Lyophilized PNGase F and *RapiGest*™ SF are included in OmniGLYZOR.

## Quick Guide

### 1 Equilibration

- Equilibrate the OmniGLYZOR Microspin column with  $3 \times 300 \mu\text{l}$  reaction buffer. Centrifuge at  $200 \times g$  for 1 min.



### 2 Deglycosylation

- Add the glycoprotein to the column and cap the column. Incubate in a thermal mixer at  $37^\circ\text{C}$  with  $650\text{-}850\text{ rpm}^1$  mixing for  $1\text{-}4\text{ h}^2$ .



### 3 Collection

- Centrifuge at  $1000 \times g$  for 1 min to collect the deglycosylated protein. For maximum recovery, add  $100 \mu\text{l}$  reaction buffer, resuspend media and centrifuge at  $1000 \times g$  for 1 min.



# PRODUCT DESCRIPTION

OmniGLYZOR contain a resin with a mixture of immobilized enzymes covalently coupled to agarose beads for removal of N- and mucin-type O-glycans on antibodies, fusion proteins and other glycosylated proteins. The enzymes are recombinantly expressed in *E. coli*.

OmniGLYZOR contain enzymes necessary for the removal of N- and simple O-linked glycans (sialyl core 1 and Tn antigen ( $\alpha$ -GalNAc)).

The following activities are included:

- PNGase F
- O-Glycosidase
- Sialidase
- $\alpha$ -GalNAcase

Removal of glycans is widely used to reduce heterogeneity to facilitate analysis of glycoproteins by for example mass spectrometry. Deglycosylation

can also be used to study the functional role of the glycans.

The glycoprotein sample is incubated with the OmniGLYZOR Microspin column for 1-4 h under native reaction conditions. The deglycosylated glycoprotein is then easily collected by a centrifugation step.

Certain N-glycosylation sites are poorly or not at all accessible to PNGaseF unless the substrate protein is denatured. Should there be any N-glycans left after incubation on the OmniGLYZOR Microspin column, they can be removed by an additional deglycosylation step under denaturing conditions using the lyophilized PNGaseF and *RapiGest* SF surfactant\* included in the kit. The enzymes hydrolyzing the O-glycans do not benefit from denaturation which is why the OmniGLYZOR Microspin columns should only be used under native conditions.

# PRODUCT DESCRIPTION

## Content and Storage

Reagents included in OmniGLYZOR:

- OmniGLYZOR Microspin columns (5x50-100 µg or 10x50-100 µg)
- Lyophilized PNGaseF (1000 units, 1 vial)
- *RapiGest* SF Surfactant (1 mg, 1 or 2 vials)

OmniGLYZOR Microspin columns contain sufficient material to remove N- and O- glycans from 50-100 µg glycoprotein. The resin is supplied in 20% EtOH with no preservatives added.

OmniGLYZOR Microspin columns are shipped cold and should be stored at +4-8°C upon arrival.

**Do not freeze the microspin columns!**

PNGaseF is supplied lyophilized in 50 mM HEPES buffer pH 7.5, with no preservatives added.

PNGaseF is shipped cold and should be stored at -20°C upon arrival. After reconstitution, PNGaseF is stable for at least 1 month at +4-8°C.

*RapiGest* SF\* should be stored at room temperature. Once reconstituted in high purity water or a buffer (pH 7-10), the solution is stable for **one week** when stored at +2-8°C.

OmniGLYZOR is for R&D use only.

*\* RapiGest™ SF Surfactant is a trademark of Waters Corporation.*

# DETAILED PROTOCOL

Use lids and bottom caps during the incubation.

Before centrifugation, remove the bottom cap and loosen the lid (do not remove the lid).

## Equipment Required

- Centrifuge for microcentrifuge tubes
- Thermal mixer compatible with microcentrifuge tubes

## Additional Materials Required

- Reaction buffer: PBS pH 7.4<sup>3</sup>
- Collection tubes: Microcentrifuge tubes (1.5-2 ml)

For preparation of denatured samples for MS analysis:

- DTT
- Formic Acid
- 8M guanidine hydrochloride pH 8.5 buffered aqueous solution, or equivalent



## Deglycosylation

### Sample Preparation

- Prepare the glycoprotein in 50-100  $\mu$ l reaction buffer per column. Recommended amount of glycoprotein is 50-100  $\mu$ g per column.

### 1 Equilibration

- Break off the bottom cap of the column (save the cap) and place the column in a collection tube. Loosen the lid.
- Centrifuge at 200  $\times$  g for 1 min to remove the storage solution.
- Equilibrate the column by adding 300  $\mu$ l reaction buffer and centrifuge at 200  $\times$  g for 1 min.
- Repeat the equilibration step an additional two times.
- Seal the spin column with the bottom cap.

### 2 Deglycosylation

- Add the glycoprotein to the column (50-100  $\mu$ g in 50-100  $\mu$ l reaction buffer).
- Seal the column with the top lid.

# DETAILED PROTOCOL

- Fully suspend the media. Do not invert the column to prevent resin getting stuck in the lid.
- Incubate the column in a thermal mixer at 37°C with 650-850 rpm<sup>1</sup> mixing for 1-4 h<sup>2</sup>.

### **3 Collection of Deglycosylated Protein**

- Remove the bottom cap and place the column in a collection tube. Loosen the top lid.
- Centrifuge the column at 1 000 × g for 1 min to recover the deglycosylated protein.
- For Maximum Recovery of the Sample:
  - Seal the spin column with the bottom cap.
  - Add 100 µl reaction buffer.
  - Seal the column and make sure the media is fully resuspended.
  - Remove the bottom cap and place the column in a collection tube. Loosen the top lid.
  - Centrifuge the column at 1 000 × g for 1 min to collect the material.
  - Pool the collected fractions.

## Additional Deglycosylation with PNGase F under Denaturing Reaction Conditions

### Sample Preparation

- Use a part or the entire collected fraction from the OmniGLYZOR column for further deglycosylation under denaturing conditions if necessary.
- Reconstitute the content of the *RapiGest* SF vial in 20 µl reaction buffer to obtain a 5 % solution. Avoid pipetting up-and-down to avoid the formation of foam.
- Add 5% *RapiGest* SF solution to the glycoprotein solution to a final *RapiGest* SF concentration of 0.1-0.2 %.
- Optional Reduction: add DTT to a final concentration of up to 50 mM to reduce disulfide bonds.
- Incubate the glycoprotein solution at 90°C for 5 min to denature the glycoprotein.
- Let the solution cool to room temperature.

# DETAILED PROTOCOL

## 1 Prepare PNGase F

- Reconstitute PNGase F in 50  $\mu$ l ddH<sub>2</sub>O to a concentration of 20 units/ $\mu$ l.

## 2 Add PNGase F

- Add 1 unit PNGase F / 1  $\mu$ g glycoprotein.

## 3 Deglycosylation

- Incubate for 30 min at 50°C<sup>2</sup>.

### **Guidelines for Preparation of Denatured Samples for MS Analysis**

#### **Acidification – hydrolysis of *RapiGest SF*:**

Add formic acid to a final concentration of 1%. Incubate at 37°C for 45 min. Slight cloudiness of the sample may be observed.

**Solubilization of the protein:** Add guanidine hydrochloride solution to a final concentration of 4 M and mix well to solubilize any precipitated protein. Centrifuge the solution at 16000  $\times$  g for 10 min. Transfer the supernatant to LC-MS vials for analysis.

## Notes

- 1. The optimal speed of mixing depends on the model of thermal mixer that is used and might be higher than the range given here. The resin must stay in suspension during the entire incubation for optimal performance.*
- 2. Incubation times depend on the glycoprotein. Some substrates might need longer incubation, especially substrates with poorly accessible N-glycans might need overnight incubation to be removed under native conditions. Alternatively, perform the additional deglycosylation step under denaturing conditions.*
- 3. Optimizations may be required if a reaction buffer other than the recommended one is used.*

## Quality Control

OmniGLYZOR is tested to meet the specifications and lot-to-lot consistency.

OmniGLYZOR is tested for the absence of microbial contamination using blood agar plates, Sabouraud dextrose agar plates and fluid thioglycollate medium.

## Related Products

### **OglyZOR®**

Hydrolysis of core 1 O-glycans

### **GalactEXO™**

Hydrolysis of  $\beta$ 1-3,4 galactose

### **GalNAcEXO™**

Hydrolysis of  $\alpha$ -linked GalNAcs

### **GlycINATOR®**

Hydrolysis of IgG Fc glycans

## Legal and Disclaimers

---

All rights reserved. Genovis products may be covered by one or more patents, trademarks and copyrights owned or controlled by Genovis AB. For more information about commercial rights, please contact the Genovis team at [licensing@genovis.com](mailto:licensing@genovis.com). Genovis products are intended for research use only. They are not intended to be used for therapeutic or diagnostic purposes in humans or animals. All goods and services are sold subject to Genovis' General Terms and Conditions of Sale.

*RapiGest*<sup>™</sup> SF Surfactant included in OmniGLYZOR. *RapiGest*<sup>™</sup> SF Surfactant is a trademark of Waters Corporation.

© Genovis AB.



## **USA & Canada**

---

Genovis Inc.  
245 First Street, Suite 1800  
Cambridge, MA 02142  
USA

Customer service: 1-617-444-8421  
Order phone (toll free): 1-855-782-0084  
Order fax: 1-858-524-3006  
Email: [orders.us@genovis.com](mailto:orders.us@genovis.com)

## **EMEA & Asia**

---

Genovis AB  
Box 790  
SE-220 07 Lund  
Sweden

Customer service: +46 46 10 12 30  
Order phone: +46 46 10 12 30  
Order fax: +46 46 12 80 20  
Email: [order@genovis.com](mailto:order@genovis.com)