

# ColdEye+




THE ONLY SCANNER FOR ALL APPLICATIONS



Operational with  
**any** sample format

Operational from  
**warm to ultra-cold**





# ColdEye+

Reliable Operation	Between +25° C and -150° C
Sample format	Compatible with any SBS Rack or Box
Barcode identification	Any 2D bottom barcode (vials) and side barcode (linear and 2D) at SBS Rack or Box
Scanning time	Approx. 3 seconds
File export	csv/xml
Cost efficiency	Integrated analysis tool to avoid additional costs (no additional software from ASKION is required)
Power supply	Integrated battery for operation up to 5h (at -150°C)
Connectivity	WLAN   LAN   HDMI   USB
Portable device	Yes
Advantages	Scanning at cold/ambient temperatures (e.g. -80/-140°C) enables best sample quality due to reduced temperature changes