




StemFit Basic04

The Next Generation Media




Animal-Origin
Free (AOF)



Weekend-Free
Feeding



Single-Cell
Expansion



One Bottle
Composition



Animal-Origin Free



Minimize risk of virus contamination and lot-to-lot variation with animal-origin free formula.



Weekend-Free Feeding




Skip feeding on weekends for flexible scheduling.

	Mon	Tue	Wed	Thu	Fri	Sat	Sun	
pattern 1	F	OFF	F	P	F	OFF	OFF	
pattern 2	F	F	P	F	F	OFF	OFF	P
pattern 3	F	P	F	OFF	F	OFF	OFF	F

P
Passage
F
Feed

Top class cost-effectiveness

Medium volume for a well





Culture cost for a well



Single-Cell Expansion



Experience superior expansion with single-cell passaging and consistent performance with various matrices.



hiPSCs (1210B2) were cultured with StemFit Basic04 on respective cell adhesion substrates. The graph shows accumulated fold expansion.

One Bottle Composition



Achieve sterilization by simplifying the procedure with ready to use packaging.



Lineup

For Basic Research

StemFit Basic02



For Clinical Research

StemFit Basic03



Frozen form, easy to use

ActivinA



Eat Well, Live Well.



For further information, please contact here. ✉ stemfit@ajinomoto.com

AJINOMOTO CO., INC. AminoScience Division <https://www.ajitrade.com/stemfit/>
15-1, Kyobashi 1-Chome, Chuo-Ku, Tokyo 104-8315, Japan