
IMCSTIPS® Catalog

When You Need Uncompromising Performance, IMCStips® Deliver

IMCStips are the definitive choice when it comes to sample preparation. Designed with efficiency and user-friendliness in mind, this product uses patented dispersive solid-phase extraction (dSPE) to automate workflows while enhancing accuracy and consistency. Unlike fixed-bed SPE devices, IMCStips contain loosely packed resins that mix with sample solutions during aspirate and dispense cycles, ensuring maximum contact between each resin and your analytes of interest.



Why Choose IMCStips?

- Consistent, high recoveries
- Customized applications
- Flexible sample volumes
- Streamlined, automated workflow

We Make Automating Sample Preparation as Easy as *Just Click Go*

You need fast, accurate results, so we designed IMCStips to make their adoption as easy as possible. In addition to on-site technical support, IMCStips come with software scripts designed for multiple automated liquid handler (ALH) platforms including Hamilton, Dynamic Devices, and INTEGRA VIAFLO 96. Through collaborations with our applications team, we create fully-developed user guidelines to help you navigate each phase of the implementation process, ensuring seamless integration of the technology into your laboratory's customized workflows.

AFFINITY HAMILTON

Protein A (MabSelect SuRe™ LX)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	5 µL	8	04T-H3R80-1-5-8
300 µL	5 µL	96	04T-H3R80-1-5-96
300 µL	5 µL	8	04T-H3R80-1A-5-8
300 µL	5 µL	96	04T-H3R80-1A-5-96
1 mL	10 µL	8	04T-H6R80-1A-10-8
1 mL	10 µL	96	04T-H6R80-1A-10-96
1 mL	25 µL	8	04T-H6R80-1A-25-8
1 mL	25 µL	96	04T-H6R80-1A-25-96
1 mL	50 µL	8	04T-H6R80-1A-50-8
1 mL	50 µL	96	04T-H6R80-1A-50-96

Protein A (MabSelect™ PrismA)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	5 µL	8	04T-H3R80A-1-5-8
300 µL	5 µL	96	04T-H3R80A-1-5-96
300 µL	5 µL	8	04T-H3R80A-1A-5-8
300 µL	5 µL	96	04T-H3R80A-1A-5-96
1 mL	25 µL	8	04T-H6R80A-1-25-8
1 mL	25 µL	96	04T-H6R80A-1-25-96
1 mL	50 µL	8	04T-H6R80A-1-50-8
1 mL	50 µL	96	04T-H6R80A-1-50-96
1 mL	100 µL	8	04T-H6R80A-1-100-8
1 mL	100 µL	96	04T-H6R80A-1-100-96

Protein A (MabSelect SuRe™ pcc)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
1 mL	25 µL	8	04T-H6R80P-1-25-8
1 mL	25 µL	96	04T-H6R80P-1-25-96
1 mL	50 µL	8	04T-H6R80P-1-50-8
1 mL	50 µL	96	04T-H6R80P-1-50-96
1 mL	100 µL	8	04T-H6R80P-1-100-8
1 mL	100 µL	96	04T-H6R80P-1-100-96

Protein A (MabSelect SuRe™)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
1 mL	10 µL	8	04T-H6R80S-1-10-8
1 mL	10 µL	96	04T-H6R80S-1-10-96
1 mL	15 µL	8	04T-H6R80S-1-15-8
1 mL	15 µL	96	04T-H6R80S-1-15-96
1 mL	25 µL	8	04T-H6R80S-1-25-8
1 mL	25 µL	96	04T-H6R80S-1-25-96
1 mL	50 µL	8	04T-H6R80S-1-50-8
1 mL	50 µL	96	04T-H6R80S-1-50-96
1 mL	75 µL	8	04T-H6R80S-1-75-8
1 mL	75 µL	96	04T-H6R80S-1-75-96

Protein A (Capture Select™ CH1-XL)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
1 mL	25 µL	8	04T-H6R83-1-25-8
1 mL	25 µL	96	04T-H6R83-1-25-96
1 mL	100 µL	8	04T-H6R83-1-100-8
1 mL	100 µL	96	04T-H6R83-1-100-96

Protein G (GammaBind™ Plus)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	5 µL	8	04T-H3R81-1-5-8
300 µL	5 µL	96	04T-H3R81-1-5-96
1 mL	25 µL	8	04T-H6R81-1A-25-8
1 mL	25 µL	96	04T-H6R81-1A-25-96
1 mL	50 µL	8	04T-H6R81-1A-50-8
1 mL	50 µL	96	04T-H6R81-1A-50-96

Protein A/G

Tip Size	Resin Amount	Tip Quantity	Catalog Number
1 mL	50 µL	8	04T-H6R82-1-50-8
1 mL	50 µL	96	04T-H6R82-1-50-96

Streptavidin (34 µm)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	4 µL	8	04T-H3R85-1-4-8
300 µL	4 µL	96	04T-H3R85-1-4-96

Streptavidin HC (45 - 105 µm)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	10 µL	8	04T-H3R85P-1-10-8
300 µL	10 µL	96	04T-H3R85P-1-10-96

Ni-IMAC (Ni Sepharose™ 6 Fast Flow)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	5 µL	8	04T-H3R72-1-5-8
300 µL	5 µL	96	04T-H3R72-1-5-96
300 µL	10 µL	8	04T-H3R72-1-10-8
300 µL	10 µL	96	04T-H3R72-1-10-96
1 mL	25 µL	8	04T-H6R72-1-25-8
1 mL	25 µL	96	04T-H6R72-1-25-96
1 mL	50 µL	8	04T-H6R72-1-50-8
1 mL	50 µL	96	04T-H6R72-1-50-96
1 mL	100 µL	8	04T-H6R72-1-100-8
1 mL	100 µL	96	04T-H6R72-1-100-96

Ni-IMAC (Ni Sepharose™ excel)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
1 mL	25 µL	8	04T-H6R72XL-1-25-8
1 mL	25 µL	96	04T-H6R72XL-1-25-96
1 mL	50 µL	8	04T-H6R72XL-1-50-8
1 mL	50 µL	96	04T-H6R72XL-1-50-96
1 mL	75 µL	8	04T-H6R72XL-1-75-8
1 mL	75 µL	96	04T-H6R72XL-1-75-96
1 mL	100 µL	8	04T-H6R72XL-1-100-8
1 mL	100 µL	96	04T-H6R72XL-1-100-96

Ni-NTA Superflow

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	10 µL	8	04T-H3R72Q-1A-10-8
300 µL	10 µL	96	04T-H3R72Q-1A-10-96
1 mL	50 µL	8	04T-H6R72Q-1-50-8
1 mL	50 µL	96	04T-H6R72Q-1-50-96

Capto Core 400

Tip Size	Resin Amount	Tip Quantity	Catalog Number
1 mL	10 µL	8	04T-H6M05-1-10-8
1 mL	10 µL	96	04T-H6M05-1-10-96

LigaG 1K (1K agarose)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	10 µL	8	04T-H3R94-1-10-8
300 µL	10 µL	96	04T-H3R94-1-10-96
1 mL	25 µL	8	04T-H6R94-1-25-8
1 mL	25 µL	96	04T-H6R94-1-25-96

LigaM 30K (30K agarose)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	10 µL	8	04T-H3R95-1-10-8
300 µL	10 µL	96	04T-H3R95-1-10-96
1 mL	25 µL	8	04T-H6R95-1-25-8
1 mL	25 µL	96	04T-H6R95-1-25-96

XT 4Flow HC (Strep-Tactin®XT 4Flow® high capacity)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	5 µL	8	04T-H3R88F-1-5-8
300 µL	5 µL	96	04T-H3R88F-1-5-96
300 µL	10 µL	8	04T-H3R88F-1-10-8
300 µL	10 µL	96	04T-H3R88F-1-10-96
300 µL	10 µL	8	04T-H3R88F-1A-10-8
300 µL	10 µL	96	04T-H3R88F-1A-10-96
1 mL	15 µL	8	04T-H6R88F-1-15-8
1 mL	15 µL	96	04T-H6R88F-1-15-96
1 mL	25 µL	8	04T-H6R88F-1-25-8
1 mL	25 µL	96	04T-H6R88F-1-25-96
1 mL	50 µL	8	04T-H6R88F-1-50-8
1 mL	50 µL	96	04T-H6R88F-1-50-96

SIZE EXCLUSION HAMILTON

SizeX$_{100}$

Tip Size	Sample Volume	Tip Quantity	Catalog Number
1 mL	100 µL	8	04T-H6R76-0A-220-8
1 mL	100 µL	96	04T-H6R76-0A-220-96

SizeX$_{150}$

Tip Size	Sample Volume	Tip Quantity	Catalog Number
1 mL	150 µL	8	04T-H6R76-0A-350-8
1 mL	150 µL	96	04T-H6R76-0A-350-96

Capto Core 400

Tip Size	Sample Volume	Tip Quantity	Catalog Number
1 mL	10 µL	8	04T-H6M05-1-10-8
1 mL	10 µL	96	04T-H6M05-1-10-96

REVERSE PHASE HAMILTON

RP

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	2 mg	8	04T-H3R05-1-2-8
300 µL	2 mg	96	04T-H3R05-1-2-96
300 µL	5 mg	8	04T-H3R05-1-5-8
300 µL	5 mg	96	04T-H3R05-1-5-96
1 mL	10 mg	8	04T-H6R05-1-10-8
1 mL	10 mg	96	04T-H6R05-1-10-96
1 mL	20 mg	8	04T-H6R05-1-20-8
1 mL	20 mg	96	04T-H6R05-1-20-96

Large Pore C$_4$

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	2 mg	8	04T-H3R41-1-2-8
300 µL	2 mg	96	04T-H3R41-1-2-96
1 mL	10 mg	8	04T-H6R41-1-10-8
1 mL	10 mg	96	04T-H6R41-1-10-96
1 mL	20 mg	8	04T-H6R41-1-20-8
1 mL	20 mg	96	04T-H6R41-1-20-96

OligoR

Tip Size	Resin Amount	Tip Quantity	Catalog Number
1 mL	50 mg	8	04T-H6R52-1-50-8
1 mL	50 mg	96	04T-H6R52-1-50-96
1 mL	50 mg	480	04T-H6R52-1-50-480

C$_8$

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	2 mg	8	04T-H3R53-1-2-8
300 µL	2 mg	96	04T-H3R53-1-2-96
1 mL	10 mg	8	04T-H6R53-1-10-8
1 mL	10 mg	96	04T-H6R53-1-10-96
1 mL	20 mg	8	04T-H6R53-1-20-8
1 mL	20 mg	96	04T-H6R53-1-20-96

C$_{18}$

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	2 mg	8	04T-H3R52-1-2-8
300 µL	2 mg	96	04T-H3R52-1-2-96
300 µL	5 mg	8	04T-H3R52-1-5-8
300 µL	5 mg	96	04T-H3R52-1-5-96
1 mL	10 mg	8	04T-H6R52-1-10-8
1 mL	10 mg	96	04T-H6R52-1-10-96
1 mL	20 mg	8	04T-H6R52-1-20-8
1 mL	20 mg	96	04T-H6R52-1-20-96

PHOSPHOPEPTIDE HAMILTON

Zirconia

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	5 mg	8	04T-H3R68-1-5-8

PolyTi

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	10 µL	8	04T-H3R73-1-10-8
300 µL	10 µL	96	04T-H3R73-1-10-96

ION EXCHANGE HAMILTON

Macro SCX

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	2 mg	8	04T-H3R30-1-2-8
300 µL	2 mg	96	04T-H3R30-1-2-96
300 µL	5 mg	8	04T-H3R30-1-5-8
300 µL	5 mg	96	04T-H3R30-1-5-96
1 mL	10 mg	8	04T-H6R30-1-10-8
1 mL	10 mg	96	04T-H6R30-1-10-96
1 mL	20 mg	8	04T-H6R30-1-20-8
1 mL	20 mg	96	04T-H6R30-1-20-96

Macro SAX

Tip Size	Resin Amount	Tip Quantity	Catalog Number
300 µL	5 mg	8	04T-H3R03R-1-5-8
300 µL	5 mg	96	04T-H3R03R-1-5-96
1 mL	20 mg	8	04T-H6R03R-1-20-8
1 mL	20 mg	96	04T-H6R03R-1-20-96

Macro WAX

Tip Size	Resin Amount	Tip Quantity	Catalog Number
1 mL	20 mg	8	04T-H6R01R-1-20-8
1 mL	20 mg	96	04T-H6R01R-1-20-96

Macro WCX

Tip Size	Resin Amount	Tip Quantity	Catalog Number
1 mL	20 mg	8	04T-H6R02R-1-20-8
1 mL	20 mg	96	04T-H6R02R-1-20-96

µSAX (Fractogel EMD TMAE)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
1 mL	25 µL	8	04T-H6R33N-1-25-8
1 mL	25 µL	96	04T-H6R33N-1-25-96

AFFINITY dynamicdevices

Ni-IMAC (Ni Sepharose™ 6 Fast Flow)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
1 mL	15 µL	8	04T-D1R72-1-15-8
1 mL	15 µL	96	04T-D1R72-1-15-96
1 mL	25 µL	8	04T-D1R72-1-25-8
1 mL	25 µL	96	04T-D1R72-1-25-96
1 mL	50 µL	8	04T-D1R72-1-50-8
1 mL	50 µL	96	04T-D1R72-1-50-96
1 mL	100 µL	8	04T-D1R72-1-100-8
1 mL	100 µL	96	04T-D1R72-1-100-96

Ni-IMAC (Ni Sepharose™ excel)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
1 mL	25 µL	8	04T-D1R72XL-1-25-8
1 mL	25 µL	96	04T-D1R72XL-1-25-96

Protein A (MabSelect SuRe™ LX)

Tip Size	Resin Amount	Tip Quantity	Catalog Number
1 mL	25 µL	8	04T-D1R80-1-25-8
1 mL	25 µL	96	04T-D1R80-1-25-96
1 mL	50 µL	8	04T-D1R80-1-50-8
1 mL	50 µL	96	04T-D1R80-1-50-96

SIZE EXCLUSION dynamicdevices

$SizeX_{100}$

Tip Size	Sample Volume	Tip Quantity	Catalog Number
1 mL	100 µL	8	04T-D1R76-0-220-8
1 mL	100 µL	96	04T-D1R76-0-220-96

$SizeX_{150}$

Tip Size	Sample Volume	Tip Quantity	Catalog Number
1 mL	150 µL	8	04T-D1R76-0-350-8
1 mL	150 µL	96	04T-D1R76-0-350-96

AFFINITY INTEGRA VIAFLO 96

Protein A (MabSelect SuRe™ LX)

Tip Size	Binding Capacity (per tip)	Tip Quantity	Catalog Number
1250 µL	1 mg IgG	8	04T-I1R80-1-25-8
1250 µL	1 mg IgG	12	04T-I1R80-1-25-12
1250 µL	1 mg IgG	96	04T-I1R80-1-25-96

Ni-IMAC (Ni Sepharose™ excel)

Tip Size	Binding Capacity (per tip)	Tip Quantity	Catalog Number
1250 µL	2 mg GFP	8	04T-I1R72XL-1-25-8
1250 µL	2 mg GFP	12	04T-I1R72XL-1-25-12


Ni-IMAC (Ni Sepharose™ 6 Fast Flow)

Tip Size	Binding Capacity (per tip)	Tip Quantity	Catalog Number
1250 µL	2 mg GFP	8	04T-I1R72-1-25-8
1250 µL	2 mg GFP	12	04T-I1R72-1-25-12
1250 µL	2 mg GFP	96	04T-I1R72-1-25-96

We provide the scripts so you can *Just Click Go!*

IMCSTIPS SEAMLESSLY INTEGRATES WITH THESE AUTOMATED LIQUID HANDLING PLATFORMS


HAMILTON



Flexible Sample Processing to Alleviate Workflow Challenges

- Processes up to 96 samples in under 30 minutes
- Selectable parameters easily configured by the user include different mixing cycles, various volumes, and other optional steps
- Combines integrated purification with Hamilton's robust technology to create fully automated methods for handling time intensive manual procedures

dynamicdevices



Simplifying Complex Purification Methods

- Save time and money by processing up to 96 samples in 90 minutes or *less*
- Greater flexibility: compatible with *both* ST and VVP pipetting tools
- Get started quickly with seamless no-code automation

INTEGRA VIAFLO 96



Affordable Robotics to Optimize Routine Pipetting Tasks

- Streamlines workflows with multiple hands-free aspirate and dispense cycles
- Consistent results, high recoveries, and increased testing efficiency
- Affordably priced and easily implemented with IMCStips into nearly any lab environment

Order IMCStips with customized resins to meet your sample processing needs!



IMCSTIPS®

We're in the Science of Reliability

IMCS addresses the growing needs of both clinical and research laboratories through advanced technologies and services designed to automate manual sample preparation, enabling scientists to focus on data analysis and experimental designs. As industry leaders, our focus is on the continuous development of novel proteomics tools to dissect protein functions and discover new biomarkers. IMCStips are utilized by top pharmaceutical companies for a wide range of applications including drug target discovery, drug target validation, drug development, drug screening, and drug manufacturing.

CONTACT US FOR *FREE* SAMPLES OF IMCSTIPS!

www.imcstips.com

inquiries@imcstips.com

+1 (888) 560-2073

IMCS, Inc
110 Centrum Drive
Irmo, SC 29063 USA

NSF-ISR
Registered to ISO 9001

Integrated Micro-Chromatography Systems Quality Assurance

IMCS
integrated micro-chromatography systems

© 2022 IMCS, Inc. All rights reserved. IMCStips® is a registered trademark of Integrated Micro-Chromatography Systems, Inc. Hamilton®, Microlab®, STAR™, and NIMBUS® are owned and/or registered by Hamilton Company in the US and/other countries. Dynamic Devices® is a registered trademark of Dynamic Devices in the United States and other countries. INTEGRA is owned by INTEGRA Biosciences AG. GammaBind™, MabSelect, MabSelect SuRe™ LX, MabSelect SuRe™ pcc, Ni Sepharose™ excel, and Sepharose™ are trademarks of GE Healthcare companies and Cytiva. Strep-Tactin®, Strep-tag®II, and related trademarks are owned by IBA Lifesciences GmbH. "CaptureSelect™ and related trademarks are owned by Thermo Fisher Scientific. All third-party trademarks are the property of their respective owners.